

PRICE: USD \$39,995

CONTACT: 1-866-330-9490 | INFO@RADIANTLASER.CA

WEBSITE: WWW.RADIANTLASER.CA

BOOK A DEMO (VIRTUAL OR IN-PERSON): (403) 461-8827

WHY CHOOSE THE RADIANT PLATINUM?

- ✓ TRIPLE-WAVELENGTH TECHNOLOGY (755 NM / 808 NM / 1064 NM) - EFFECTIVELY TREATS A FULL SPECTRUM OF HAIR TYPES AND SKIN TONES.
- ✓ SAFE FOR ALL SKIN TYPES (FITZPATRICK I-VI) INCLUDING TANNED AND DARKER COMPLEXIONS.
- ✓ FAST, COMFORTABLE TREATMENTS WITH ADVANCED COOLING + ERGONOMIC HANDPIECE - IDEAL FOR HIGH-VOLUME CLINICS.
- ✓ BROAD TREATMENT RANGE - FROM DELICATE FACIAL ZONES TO LEGS, BACK, CHEST, BIKINI LINE.
- ✓ INGROWN HAIR REDUCTION - ADDRESSES AND MINIMIZES INGROWN HAIRS FOR SMOOTHER SKIN.
- ✓ SKIN REJUVENATION BENEFITS - MILD COLLAGEN STIMULATION AND IMPROVEMENT OF SKIN TONE OVER REPEATED SESSIONS.

TECHNICAL SPECIFICATIONS & FEATURES

WAVELENGTHS: 755 NM, 808 NM, 1064 NM

SCREEN SIZE: 14 INCH 4K TRUE COLOR TOUCH SCREEN

PULSE DURATIONS: 1 - 300MS

FREQUENCY: 1-10HZ

SPOT SIZE: 12*24MM OR 12*36MM

COOLING: WATER COOLING + AIR COOLING + COMPRESSOR

VOLTAGE: AC220V/110V 50HZ/60HZ

DIMENSIONS / WEIGHT: 125*56*71CM

DISPLAY / INTERFACE: TOUCHSCREEN

COMPATIBLE SYSTEM: COMPATIBLE WITH APPLE AND ANDROID PHONES

HEALTH CANADA APPROVED, SAFETY & REGULATORY COMPLIANCE

ONE KEY START, FAST OPERATION

LOWEST TEMPERATURE COULD REACH -35C

RAPID COOLING, ONLY 30S

REALLY PAINLESS WHEN MAKE THE TREATMENT

PROTECT THE CUSTOMERS SKIN

FAC TECHNOLOGY, LASER REACHS TOOT OF HAIR

ICE PAINLESS TREATMENT FEELING

COMFORTABLE EXPERIENCE

CONTINUE COOLING

ADJUSTABLE COOLING LEVEL

REAL CLINIC ADVANTAGES

HIGH THROUGHPUT

TREAT MORE CLIENTS PER DAY WITH SPEED AND COMFORT.

VERSATILITY

ONE MACHINE HANDLES MULTIPLE SKIN/HAIR TYPES.

RELIABILITY & SUPPORT

RADIANT BACKS ITS DEVICES WITH SERVICE, PARTS, AND WARRANTY SUPPORT.

REVENUE GENERATOR

ATTRACT NEW CLIENTS, UPSELL PACKAGES, AND GROW YOUR MEDSPA/CLINIC BUSINESS.



Health
Canada



EXCLUSIVE PRICING, BONUS SERVICE & SUPPLY BUNDLES

AVAILABLE ONLY FOR SHOW ATTENDEES. TALK TO OUR REPS ON THE FLOOR TO UNLOCK YOUR SPECIAL PACKAGE.

FINANCING & PURCHASE OPTIONS

WE UNDERSTAND ACQUIRING CAPITAL EQUIPMENT CAN BE A HURDLE - THAT'S WHY RADIANT OFFERS FLEXIBLE FINANCING SOLUTIONS TO GET YOU UP AND RUNNING. WHETHER YOU PREFER LEASE, PAYMENT PLANS, OR VENDOR FINANCING, WE'LL HELP YOU STRUCTURE A PLAN THAT FITS YOUR BUSINESS CASH FLOW.

ASK YOUR SALES REP OR CONTACT US AT 1-866-330-9490 FOR DETAILS ON FINANCING TERMS, RATES, AND APPROVAL PROCESS.

BOOK A DEMO – SEE IT IN ACTION

VIRTUAL DEMO:

WATCH A LIVE DEMONSTRATION OF THE RADIANT PLATINUM FROM YOUR CLINIC OR OFFICE VIA VIDEO CONFERENCE. SPEAK WITH OUR LASER EXPERTS, ASK QUESTIONS, AND VIEW PERFORMANCE IN REAL TIME. CALL 403-461-8827 TO SCHEDULE.

IN-PERSON DEMO / TRAINING:

WATCH A LIVE DEMONSTRATION OF THE RADIANT PLATINUM FROM YOUR CLINIC OR OFFICE VIA VIDEO CONFERENCE. SPEAK WITH OUR LASER EXPERTS, ASK QUESTIONS, AND VIEW PERFORMANCE IN REAL TIME. CALL 403-461-8827 TO SCHEDULE.

WE CAN BRING A UNIT TO YOUR LOCATION (WITHIN SERVICE REGIONS) OR HOST YOU AT ONE OF OUR SHOWROOMS. EXPERIENCE FULL OPERATION, PERFORM SAMPLE TREATMENTS, AND RECEIVE HANDS-ON TRAINING GUIDANCE. CALL 403-461-8827 TO ARRANGE YOUR ON-SITE OR IN-CLINIC DEMO.

ADDITIONAL INFO & LEGAL/WARRANTY NOTES

- ALL DEVICE WARRANTIES, SERVICE CONTRACTS, AND MAINTENANCE OPTIONS WILL BE DETAILED IN YOUR PURCHASE AGREEMENT.
- RADIANT'S SUPPORT AND PARTS AVAILABILITY ENSURE YOU REMAIN OPERATIONAL WITH MINIMAL DOWNTIME.
- THE RADIANT PLATINUM IS HEALTH CANADA APPROVED FOR SAFETY AND PERFORMANCE ACROSS SKIN TYPES.
- RADIANT ALSO OFFERS TRADE-IN AND UPGRADE PROGRAMS; ASK YOUR REP FOR ELIGIBILITY.

HOW TO ORDER / NEXT STEPS

CONTACT SALES: AT 1-866-330-9490 OR INFO@RADIANTLASER.CA

REQUEST DEMO

SECURE FINANCING / SIGN AGREEMENT

SCHEDULE DELIVERY, INSTALLATION & TRAINING

LAUNCH YOUR SERVICES & PROMOTE YOUR NEW PREMIUM LASER OFFERING

**DON'T MISS YOUR CHANC. ATTEND A DEMO,
EXPLORE TRADE SHOW SPECIALS, AND BEGIN ELEVATING YOUR
CLINIC'S LASER HAIR REMOVAL OFFERING.**

CONTACT US TODAY:

RADIANT TECHNOLOGY CORP.

1-866-330-9490 | INFO@RADIANTLASER.CA | WWW.RADIANTLASER.CA

DEMO & TRAINING BOOKING: (403) 461-8827

