

PRICE: \$55,000 CAD

BOOK A DEMO (VIRTUAL OR IN-PERSON): (403) 461-8827

CONTACT: 1-866-330-9490 | INFO@RADIANTLASER.CA | WWW.RADIANTLASER.CA

## WHY RADIANT GENESIS STANDS OUT

- ✓ PAINLESS TECHNOLOGY + ADVANCED COOLING: DESIGNED TO DELIVER EFFECTIVE TREATMENTS WITH MINIMAL DISCOMFORT.
- ✓ SAFE FOR ALL SKIN TYPES & TONES: VERSATILE ACROSS FITZPATRICK I-VI SKIN.
- ✓ MULTI-FUNCTION PLATFORM: HAIR REMOVAL, SKIN REJUVENATION, PIGMENT CORRECTION, VASCULAR THERAPY, ACNE, STRETCH MARK REDUCTION, ROSACEA, BROKEN CAPILLARIES, INGROWN HAIR TREATMENT.
- ✓ FASTER TREATMENTS, MINIMAL DOWNTIME: OPTIMIZED FOR HIGH THROUGHPUT IN CLINIC SETTINGS.
- ✓ CLINICALLY PROVEN & RELIABLE: UTILIZES A SPECTRUM OF 530-1200 NM AND 640-1200 NM WAVELENGTHS TO BALANCE DEPTH AND ABSORPTION.

## TECHNICAL SPECIFICATIONS & KEY FEATURES

WAVELENGTH RANGES	530 - 1200 NM (SR HANDLE), 640 - 1200 NM (HR HANDLE)
SPOT SIZE	10 × 50 MM
ENERGY DENSITY	3 - 25 J/CM <sup>2</sup> (1 J STEPS)
PULSE WIDTH	1 - 20 MS IN 1 MS STEPS
PULSE INTERVAL (SUB)	10 MS (FIXED)
PULSE STRING INTERVAL	1 - 4 S
HANDLES	2 HANDLES (HR & SR)
COOLING METHOD	AIR / WATER / SEMICONDUCTOR COOLING, HANDLE COOLING TEMP 0 - 15 °C
CONTROLLER / INTERFACE	15-INCH TRUE-COLOR TOUCHSCREEN
PACKAGING DIMENSIONS	MAIN: 500 × 471 × 1085 MM; PACKAGED: 600 × 610 × 1170 MM
NET WEIGHT	MAIN UNIT: - 73 KG; PACKAGED: -95 KG
INPUT POWER	2500 VA

## ADVANTAGES IN CLINIC USE

- ✓ VERSATILITY & REVENUE EXPANSION  
OFFERING MULTIPLE TREATMENT TYPES FROM ONE DEVICE ATTRACTS BROADER CLIENTELE.
- ✓ HIGH EFFICIENCY  
LARGE SPOT SIZE AND EFFICIENT COOLING ALLOW YOU TO TREAT MORE CLIENTS PER SHIFT.
- ✓ PATIENT COMFORT & RETENTION  
THE NEAR-PAINLESS EXPERIENCE ENCOURAGES RETURN VISITS AND POSITIVE REFERRALS.
- ✓ SUPPORT & WARRANTY  
RADIANT BACKS THE GENESIS WITH SERVICE, PARTS, AND TRAINING.



**UNLOCK EXCLUSIVE TRADE SHOW PRICING, BONUS CONSUMABLES OR SERVICE PACKAGES, AND DEFERRED PAYMENT OPTIONS**

AVAILABLE ONLY DURING THE SHOW OR INITIAL ORDER WINDOW. SPEAK WITH OUR REPS TO CLAIM YOUR SPECIAL PACKAGE.

**FLEXIBLE FINANCING OPTIONS**

WE UNDERSTAND CAPITAL EQUIPMENT IS A SIGNIFICANT INVESTMENT. RADIANT OFFERS CUSTOM FINANCING, LEASE ARRANGEMENTS, OR PAYMENT PLANS TO HELP YOU INTEGRATE THE GENESIS INTO YOUR PRACTICE WITHOUT STRAINING CASH FLOW. **CALL 1-866-330-9490 TO EXPLORE TAILORED PAYMENT TERMS AND FINANCING APPLICATIONS.**

**TRY BEFORE YOU BUY – BOOK YOUR DEMO**

**VIRTUAL DEMO:** SEE THE GENESIS IN ACTION LIVE VIA VIDEO. ASK QUESTIONS, OBSERVE TREATMENTS IN REAL TIME, AND GET A FULL WALKTHROUGH. CALL 403-461-8827 TO SCHEDULE.

**IN-PERSON / ONSITE DEMO & TRAINING:**

WE CAN BRING A UNIT TO YOUR CLINIC (WHERE SERVICEABLE) OR HOST YOU IN ONE OF OUR FACILITIES. HANDS-ON OPERATION, TEST TREATMENTS, AND TRAINING INCLUDED. CALL 403-461-8827 TO SCHEDULE.

WE CAN BRING A UNIT TO YOUR LOCATION (WITHIN SERVICE REGIONS) OR HOST YOU AT ONE OF OUR SHOWROOMS. EXPERIENCE FULL OPERATION, PERFORM SAMPLE TREATMENTS, AND RECEIVE HANDS-ON TRAINING GUIDANCE. CALL 403-461-8827 TO ARRANGE YOUR ON-SITE OR IN-CLINIC DEMO.

**ADDITIONAL DISCLOSURES & SUPPORT**

- ALL WARRANTIES, SERVICE CONTRACTS, AND MAINTENANCE TERMS ARE DETAILED IN THE PURCHASE AGREEMENT.
- RADIANT PROVIDES SUPPORT, PARTS, AND SERVICE TO MINIMIZE DOWNTIME.
- RADIANT ALSO OFFERS TRADE-IN, UPGRADE, AND BUNDLE PROGRAMS - INQUIRE WITH YOUR SALES REP.
- THE RADIANT GENESIS IS CERTIFIED UNDER HEALTH CANADA AND FDA FOR SAFETY AND CLINICAL PERFORMANCE.

**HOW TO ORDER / NEXT STEPS**

CONTACT SALES: AT 1-866-330-9490 OR [INFO@RADIANTLASER.CA](mailto:INFO@RADIANTLASER.CA)

REQUEST DEMO: VIRTUAL OR IN PERSON AT 403-461-8827

FINALIZE FINANCING & AGREEMENT

SCHEDULE DELIVERY, INSTALLATION & TRAINING

LAUNCH YOUR SERVICES & PROMOTE

**DON'T WAIT – ELEVATE YOUR CLINIC'S OFFERING WITH THE GENESIS. BOOK YOUR DEMO NOW, SECURE SHOW SPECIALS, AND AMPLIFY YOUR PATIENT OUTCOMES.**

**CONTACT US TODAY:**

RADIANT TECHNOLOGY CORP.

1-866-330-9490 | [INFO@RADIANTLASER.CA](mailto:INFO@RADIANTLASER.CA)

DEMO & TRAINING BOOKING: (403) 461-8827

WEBSITE: [WWW.RADIANTLASER.CA](http://WWW.RADIANTLASER.CA)

