

PRICE: \$24,000 CAD

BOOK A DEMO (VIRTUAL OR IN-PERSON): 1-866-330-9490

EMAIL: [INFO@RADIANTLASER.CA](mailto:INFO@RADIANTLASER.CA) | WEBSITE: [WWW.RADIANTLASER.CA](http://WWW.RADIANTLASER.CA)

### PROFESSIONAL FACIAL IMAGING + 10-MODE SKIN PERFORMANCE ANALYSIS

Meet Radiant Lumina™ – the next evolution in intelligent skin diagnostics. Designed for high-performance clinics and luxury aesthetic studios, Lumina combines advanced multi-spectral imaging with intuitive software to reveal what the eye can't see, turning every consultation into a precise, confidence-building roadmap for transformation.

From sun damage hiding beneath the surface to pore congestion and redness patterns, Radiant Lumina delivers a 360° skin story in seconds, helping you communicate results clearly, recommend treatments accurately, and track change over time with undeniable visual proof.

### WHAT RADIANT LUMINA DOES?

Radiant Lumina captures high-resolution facial images using standard, cross-polarized, and UV-style illumination to analyze skin from surface to sub-surface layers. The system then generates an easy-to-read report with visual overlays, comparative scoring, and progress tracking.

Think of it as your clinic's truth mirror – objective, consistent, and impressively detailed.

### KEY FEATURES

- ✓ MULTI-SPECTRAL FACIAL IMAGING (SURFACE + SUBSURFACE VISIBILITY)
- ✓ HIGH-RESOLUTION CAPTURE FOR FULL-FACE AND ZONE ANALYSIS
- ✓ FAST, GUIDED SOFTWARE WORKFLOW
- ✓ CLIENT-FRIENDLY REPORTS WITH OVERLAYS & SCORES
- ✓ TREATMENT TRACKING DASHBOARD
- ✓ ELEGANT, COMPACT HARDWARE DESIGN
- ✓ IDEAL FOR MED-SPAS, DERMATOLOGY, PMU STUDIOS, AND COSMETIC CLINICS



## 10 MODES OF SKIN DIAGNOSTIC PERFORMANCE

Radiant Lumina includes 10 powerful analysis modes, each highlighting a key aspect of skin health and aging:

### Wrinkle Mapping

Detects fine lines and deeper wrinkle pathways, showing depth and pattern progression.

### Texture & Smoothness Index

Evaluates micro-roughness, unevenness, and skin grain for refinement planning.

### Pore Clarity & Size Analysis

Measures visible pores and congestion patterns to guide extraction, resurfacing, and tightening protocols.

### UV / Hidden Sun Damage Scan

Reveals underlying photodamage not visible in normal light - ideal for prevention education and sun-repair plans.

### Brown Spot & Pigment Distribution

Identifies pigment clusters, age spots, and uneven melanin activity across facial zones.

### Redness & Vascular Pattern Detection

Maps diffuse redness and capillary-based inflammation for rosacea, sensitivity, or post-procedure monitoring.

### Porphyryn / Bacterial Activity Reading

Highlights bacterial buildup in pores associated with acne and breakouts.

### Spot & Blemish Clarity Tracking

Quantifies visible spots and active blemish zones for corrective treatment and skincare planning.

### Skin Age Performance Score

Provides an estimated “skin age” performance metric based on the above categories - powerful for motivation and retention.

### Before-&-After Progress Comparison

Tracks changes across visits with side-by-side imaging, percent improvements, and treatment timeline overlays.

## CONTACT US TODAY:

RADIANT TECHNOLOGY CORP.

PHONE: 1-866-330-9490

EMAIL: [INFO@RADIANTLASER.CA](mailto:INFO@RADIANTLASER.CA)

WEBSITE: [WWW.RADIANTLASER.CA](http://WWW.RADIANTLASER.CA)

