

PRICE: \$69,000 CAD

BOOK A DEMO (VIRTUAL OR IN-PERSON): 1-866-330-9490

EMAIL: INFO@RADIANTLASER.CA | WEBSITE: WWW.RADIANTLASER.CA

THE FUTURE OF SKIN TRANSFORMATION HAS ARRIVED

Experience the cutting-edge precision of the Radiant Quantum Picosecond Laser - a next generation device engineered for unparalleled performance in advanced skin rejuvenation, pigmentation removal, and tattoo clearance. Using ultra-short picosecond pulse technology, this laser delivers energy to the skin in trillionths of a second, shattering pigment particles without damaging surrounding tissue.

FASTER RESULTS, MINIMAL DOWNTIME, AND MAXIMUM PATIENT COMFORT.

The Radiant Quantum is designed for professionals who demand both power and precision. With multiple wavelengths and customizable settings, it effectively treats a wide range of indications including:

- ✓ SUN DAMAGE, AGE SPOTS, AND MELASMA
- ✓ ACNE SCARS AND FINE LINES
- ✓ MULTI-COLORED TATTOO REMOVAL
- ✓ SKIN TIGHTENING AND COLLAGEN REMODELING

WHERE LIGHT MEETS PERFECTION.






Featuring advanced cooling technology and intelligent energy calibration, this system ensures consistent, safe, and effective treatments for all skin types.

SPECIFICATION

POWER	2500W
LASER TYPE	1064NM FOCUS LENS PROBE
FREQUENCY	2-10HZ
DURATION	DYNAMIC LIGHT 750PS
ENERGY	1064@100MJ-1000MJ,STEP 200MJ, 532@50MJ-500MJ, STEP 100MJ
ENERGY STABILITY	3%@1064NM/532NM (STANDARD DEVIATION/AVERAGE)
COOLING WAY	AIR COOLING + WATER COOLING
DIMENSION/WEIGHT	MACHINE 98*36*88CM/100KG PACKAGE 100*47*113/141KG



ALL PIGMENTATION CONCERNS



HOW TO ORDER / NEXT STEPS

CONTACT SALES / REQUEST DEMO: AT 1-866-330-9490 OR INFO@RADIANTLASER.CA

FINALIZE FINANCING & AGREEMENT

SCHEDULE DELIVERY, INSTALLATION & TRAINING

LAUNCH YOUR SERVICES & PROMOTE

CONTACT US TODAY:

RADIANT TECHNOLOGY CORP.

1-866-330-9490 | INFO@RADIANTLASER.CA

WEBSITE: WWW.RADIANTLASER.CA

CLINIC DATA

A course of treatment is 3-10 times, according to pigmentation's age/content/pigment type.

THE REACTION OF PIGMENT UNDER SKIN MIRROR



BEFORE

AFTER 1ST
TREATMENT

AFTER 4TH
TREATMENT

AFTER 5TH
TREATMENT

MALE / AGE 27 /
III SKIN TYPE

PIGMENT TYPE TATTOO
TATTOO AGE 5 YEARS
TATTOO COLOR BLUE AND GREEN
TREATMENT TIMES 5 TIMES

FEMALE / AGE 41 /
III SKIN TYPE

PIGMENT TYPE CHLOASMA
PIGMENT AGE 5 YEARS
PIGMENT COLOR BROWN
TREATMENT TIMES 6 TIMES



FEMALE / AGE 42 /
II SKIN TYPE

APPLICATION FINE WRINKLE REMOVAL
WRINKLE AGE 2 YEARS
TREATMENT TIMES 2 TIMES

PHOTO GALLERY

BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER



BEFORE

AFTER

